

WATER DETERMINATION IN CHALLENGING SAMPLES

NEW! FREE OF ALCOHOLS (FA) HYDRANAL REAGENTS



Hydranal™ NEXTGEN Coulomat A-FA & C-FA are the first alcohol-free reagents for Karl Fischer (KF) titration that suppress side effects and provide consistent results for many challenging alcohol-sensitive samples.

Many samples, like lithium-ion battery (LiB) electrolytes or those with additives, undergo side reactions with alcohols in common KF reagents. This leads to an increase in drift, and delayed or no titration endpoint, rendering incorrect water content results. It can even make water determination impossible.

New Hydranal alcohol-free reagents suppress these side effects and allow accurate water determination, even in difficult samples. Additionally, the new formulation does not contain CMR (carcinogenic, mutagenic and reprotoxic) substances or halogenated hydrocarbons.

NEW HYDRANAL FA-TYPE REAGENTS

- Free of alcohols
- Recommended for Li-ion battery electrolytes
- No CMR classification
- High performance
- High accuracy

Honeywell

Lithium-ion battery electrolytes are sensitive to humidity. Explosion of cell phone batteries is a known example of this. However, it is not easy to perform water determination in these samples.

Alcohol-based reagents for KF titration, like Hydranal-Coulomat AK or competitive Reagent M, are often unable to titrate LiB electrolyte samples.

Common electrolyte additives, such as vinylene carbonate (VC), fluoroethylene carbonate (FEC) or borates (e.g. lithium bis(oxalate)borate, LiBOB), cause strong side reactions leading to increased drift, missing end points and inaccurate results.

Hydranal NEXTGEN Coulomat A-FA & C-FA suppress side reactions and provide consistent performance and high accuracy even for difficult samples (see comparison on Fig. 1 and 2).

SAMPLE	COULOMAT A-FA + C-FA	COULOMAT AK + CG-K	REAGENT M
LiB electrolyte with VC/FEC	+	o	o
LiB electrolyte with LiBOB	+	-	-
LiB electrolyte with LiDFOB	+	-	-
Pure VC	+	-	-
Pure FEC	+	o	o
Pure LiBOB	+	-	-
Pure LiDFOB	+	-	-

- + accurate titration, suppressed side effects
- o limited titration (possible only for small amount of sample)
- no accurate titration possible (side reactions, increasing drift or inaccurate results)

VC - vinylene carbonate; FEC - fluoroethylene carbonate; LiBOB - lithium bis(oxalate)borate; LiDFOB - lithium difluoro(oxalato)borate

PRODUCT NUMBER	PRODUCT NAME	DESCRIPTION	PACKAGING
34471	HYDRANAL NEXTGEN Coulomat A-FA	Analyte for coulometric Karl Fischer titration of ketones and Li-ion battery electrolytes, preferred for cells with diaphragm, acetonitrile-based, free of alcohols (free of CMR substances)	500 mL
34470	HYDRANAL NEXTGEN Coulomat C-FA	Catholyte for coulometric Karl Fischer titration of ketones and Li-ion battery electrolytes, acetonitrile-based, free of alcohols (free of CMR substances)	10 x 5 mL ampoules

Although Honeywell International Inc. believes that the information contained herein is accurate and reliable, it is presented without guarantee or responsibility of any kind and does not constitute any representation or warranty of Honeywell International Inc., either expressed or implied. A number of factors may affect the performance of any products used in conjunction with user's materials, such as other raw materials, application, formulation, environmental factors and manufacturing conditions among others, all of which must be taken into account by the user in producing or using the products. The user should not assume that all necessary data for the proper evaluation of these products are contained herein. Information provided herein does not relieve the user from the responsibility of carrying out its own tests and experiments, and the user assumes all risks and liabilities (including, but not limited to, risks relating to results, patent infringement, regulatory compliance and health, safety and environment) related to the use of the products and/or information contained herein.

For more information visit:

hydranal-honeywell.com

Hydranal Center of Excellence

Tel: +49 5137 999 353
 Fax: +49 5137 999 698
 hydranal@honeywell.com



FL-21-02-ENG | 03/21
 © 2021 Honeywell International Inc.




Fig. 1. Titration of 1 M LiPF₆ in EC/DMC with 5 wt% LiBOB.




Fig. 2. Titration of pure VC.

