

# DWK Life Sciences WHEATON® COMPLETEPAK™



DESIGNED BY YOU.  
CRAFTED BY US.



DURAN  
WHEATON  
KIMBLE

Excellence in your hands

 **fisher scientific**  
part of Thermo Fisher Scientific

# WHEATON® COMPLETEPAK

The WHEATON® COMPLETEPAK product line is a range of crimp top vials, rubber stoppers and seals that are ready-to-use and are off-the-shelf, which can be customized for any application. This reduces our customer's supply chain due to our one source solution capability. Each kit comes with specific United States Pharmacopeia (USP) certificates showing that the product meets or exceeds critical USP standards that are enforced by the Federal Drug Administration (FDA). This kind of solution allows for our customers to focus on their core competency which is manufacturing, compounding and developing drug products.

## Target Applications

- Compounding pharmaceuticals
- Aseptic filling of injectable drug products
- Drug discovery
- Drug research
- Freeze drying / lyophilized drug products
- Large and small molecule drug packaging

## Key Benefits

- Sterilized, off-the-shelf, ready-to-use (RtU) products
- Variety of Kitted and stand-alone sterile components — a solution for any application
- Kitted solutions feature equal amounts of RtU components for reduced scrap rates
- Complete manufacturing & certificate of analysis traceability
- USP <788> Particulate Matter in Injections, USP <85> The Bacterial Endotoxins Test and USP <71> The Sterility Test
- CompletePAK Amber for light sensitive content packaging
- CompletePAK ModPAK Bulk configurations
- 18 month shelf life from date of processing
- Container Closure Integrity (CCI) certification available upon request

## Materials

- Vials are designed for parenteral applications and manufactured from low extractable borosilicate glass
- All stoppers are manufactured from an ultrapure bromobutyl latex free formulation with extremely high chemical purity and low gas permeability
- Omniflex stoppers offer low levels of residual binding to maximize extraction of active pharmaceutical ingredients
- Omniflex stopper elasticity delivers a more effective barrier to optimize container closure integrity and content protection
- Igloo lyophilized stoppers are suited for biologics and content with sensitivities to silicon



## ENHANCING THE EFFICIENCY OF DRUG DISCOVERY AND DRUG COMPOUNDING



## TOP REASONS WHY YOU SHOULD CHOOSE COMPLETEPAK FOR YOUR APPLICATION:

1. Products supplied ready-to-use to improve workflow
2. Single source options reduce lead times & Simplify administrative effort
3. Focuses your efforts on critical core competencies
4. Eliminates capital investments in sterile processing equipment
5. Frees up valuable clean room space
6. Validated processes & Certificate of Analysis Traceability assures quality production outcomes
7. Saves valuable time and money while maintaining a high quality product and streamlining process

## COMPLETEPAK Kits

Clear Kit Cat. No.	Amber Kit Cat. No.	Description (pieces/kit)	Qty / Case	ModPAK Cat. No. (Clear Only)	Cases / ModPAK
06-451-220	15-112-150	2mL Sterile Vial (210), 13mm Ultra Pure Sterile Serum Stopper (230), Sterile Blue Seal (230)	1	15-112-132	6
06-451-221	15-112-151	2mL Sterile Vial (210), 13mm Omni Flex 3G Sterile Serum Stopper (230), Sterile Blue Seal (230)	1	15-112-133	6
06-451-222	15-112-152	2mL Sterile Vial (210), 13mm OmniFlex 3G Sterile Igloo Lyo Stopper (230), Sterile Blue Seal (230)	1	15-112-134	6
06-451-223	15-112-153	2mL Sterile Vial (210), 13mm Ultra Pure Sterile Serum Stopper (230), Sterile Red Seal (230)	1	15-112-135	6
06-451-224	15-112-154	2mL Sterile Vial (210), 13mm Omni Flex 3G Sterile Serum Stopper (230), Sterile Red Seal (230)	1	15-112-136	6
06-451-225	15-112-155	2mL Sterile Vial (210), 13mm OmniFlex 3G Sterile Igloo Lyo Stopper (230), Sterile Red Seal (230)	1	15-112-137	6
06-451-226	15-112-156	5mL Sterile Vial (200), 20mm Ultra Pure Sterile Serum Stopper (220), Sterile Blue Seal (220)	1	15-112-138	6
06-451-227	15-112-157	5mL Sterile Vial (200), 20mm OmniFlex 3G Sterile Serum Stopper (220), Sterile Blue Seal (220)	1	15-112-139	6
06-451-228	15-112-158	5mL Sterile Vial (200), 20mm OmniFlex 3G Igloo Lyo Stopper (220), Sterile Blue Seal (220)	1	15-112-140	6
06-451-229	15-112-159	5mL Sterile Vial (200), 20mm Ultra Pure Sterile Serum Stopper (220), Sterile Red Seal (220)	1	15-112-141	6
06-451-230	15-112-160	5mL Sterile Vial (200), 20mm OmniFlex 3G Sterile Serum Stopper (220), Sterile Red Seal (220)	1	15-112-142	6
06-451-231	15-112-161	5mL Sterile Vial (200), 20mm OmniFlex 3G Sterile Igloo Lyo Stopper (220), Sterile Red Seal (220)	1	15-112-143	6
06-451-232	15-112-162	10mL Sterile Vial (158), 20mm Ultra Pure Sterile Serum Stopper (220), Sterile Blue Seal (220)	1	15-112-144	6
06-451-233	15-112-163	10mL Sterile Vial (158), 20mm OmniFlex 3G Sterile Serum Stopper (220), Sterile Blue Seal (220)	1	15-112-145	6
06-451-234	15-112-164	10mL Sterile Vial (158), 20mm OmniFlex 3G Sterile Igloo Lyo Stopper (220), Sterile Blue Seal (220)	1	15-112-146	6
06-451-235	15-112-165	10mL Sterile Vial (158), 20mm Ultra Pure Sterile Serum Stopper (220), Sterile Red Seal (220)	1	15-112-147	6
06-451-236	15-112-166	10mL Sterile Vial (158), 20mm OmniFlex 3G Sterile Serum Stopper (220), Sterile Red Seal (220)	1	15-112-148	6
06-451-237	15-112-167	10mL Sterile Vial (158), 20mm OmniFlex 3G Sterile Igloo Lyo Stopper (220), Sterile Red Seal (220)	1	15-112-149	6

## COMPLETEPAK Sterile Components\*

Cat. No.	Description	Qty / Case
<b>Clear Vials</b>		
06-451-238	2mL Clear Sterile Serum Vial	210
06-451-239	5mL Clear Sterile Serum Vial	200
06-451-240	10mL Clear Sterile Serum Vial	158
<b>Amber Vials</b>		
15-112-168	2mL Amber Sterile Serum Vial	210
15-112-169	5mL Amber Sterile Serum Vial	200
15-112-170	10mL Amber Sterile Serum Vial	158

Cat. No.	Description	Qty / Case
<b>Stoppers</b>		
06-451-241	13mm OmniFlex 3G Sterile Serum Stopper	230
06-451-242	13mm Ultra Pure Sterile Serum Stopper	230
06-451-243	13mm OmniFlex 3G Sterile Igloo Lyo Stopper	230
06-451-244	20mm OmniFlex 3G Sterile Serum Stopper	220
06-451-245	20mm Ultra Pure Sterile Serum Stopper	220
06-451-246	20mm OmniFlex 3G Sterile Igloo Lyo Stopper	220
<b>Seals</b>		
06-451-247	13mm Sterile Flip Red Seal	230
06-451-248	13mm Sterile Flip Blue Seal	230
06-451-249	20mm Sterile Flip Red Seal	220
06-451-250	20mm Sterile Flip Blue Seal	220



WHEATON® is a registered trademark of DWK Life Sciences LLC.  
Copyright © 2018 DWK Life Sciences LLC. FBRO\_018 12/2019

**In the United States**

Order online: [fishersci.com](http://fishersci.com)

Fax an order: 1-800-926-1166

Call customer service: 1-800-766-7000

**In Canada**

Order online: [fishersci.ca](http://fishersci.ca)

Fax an order: 1-800-463-2996

Call customer service: 1-800-234-7437



© 2019 Thermo Fisher Scientific Inc. All rights reserved.  
Trademarks used are owned as indicated at [fishersci.com/trademarks](http://fishersci.com/trademarks).

BN20171340