



Fisherbrand
QUALITY. RELIABILITY. VALUE.

Centrifugation

Mini and Midi Centrifuges

Mini Centrifuge

The Fisherbrand Mini Centrifuge is ideal for quick spin downs, micro-filter cell separations, PCR* and HPLC samples. It is compact and easy to use and comes as standard with a host of accessories



Key features & benefits:

- Small footprint and quiet operation
- UV resistant housing
- SnapSpin™ rotors for quick and easy change
- Compact ergonomic design and footprint
- Max speed 6,000rpm/2,000xg
- Weighted rubberised base eliminates potential movement and further dampens any noise output
- One year warranty

Technical Specification

Capacity rotor 1	6 x 1.5/2mL tubes
Capacity rotor 2	2 x 8 0.2mL strips or 16 single tubes
Adapters for rotor 1	6 x 0.5mL and 6 x 0.2mL tubes
Speed range, rpm;	6,000
RCF:	2,000xg
Dimensions:	153 x 128 x 104mm
Electrical requirements:	220-230V, 50/60Hz

Cat. No

Description

15358266

Mini centrifuge, includes rotors, tube rack, storage case, O-rings, tube adapters and low voltage, double insulated power adapter with four interchangeable plugs

For the full range of Fisherbrand products, visit www.eu.fishersci.com/go/fisherbrand

*PCR is a registered trademark of Hoffmann LaRoche Inc



Midi Centrifuge

Achieve maximal speed with minimum space. The Fisherbrand Midi Centrifuge is ideal for protocols calling for fast spins with a higher RCF requirement. Ideal for micro-filter cell separations, PCR*, clinical applications and HPLC protocols

Key features & benefits:

- Compact footprint and brush-less motor for quiet operation and low vibration
- Select and set the digital speed display in RCF or RPM
- Sample imbalance and tilt safety cut-out
- Safety interlock lid prevents opening while rotor is turning
- Back lit display showing time, speed and operational information on the same screen
- Quick spin button for short spin requirements
- UV resistant plastic, increased useful life under harsh light conditions
- Push-button lid release
- Manual override allows user to open lid in case of power failure
- Two year warranty

Supplied with two rotors:

- 12 place standard rotor with lid for 1.5/2.0mL tubes
- Eight place strip rotor for four 0.2mL tube strips or 32 single tubes



Technical Specification

Speed range, rpm;	500 to 12,500
RCF:	15 to 9,800xg
Dimensions:	203 x 171 x 114mm
Electrical requirements:	100-240V, 50/60Hz

Cat. No.	Description
----------	-------------

12972041	Midi centrifuge, includes rotors, 12 x 0.2mL and 12 x 0.5mL tube adapters, universal double insulated power adapter with four interchangeable plugs
-----------------	---

Terms and Conditions Please note that all sales are subject to our standard terms and conditions. Full details can be found on the inside back cover of our catalogues or visit our website www.fisher.co.uk

© 2015 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at www.fishersci.com/trademarks.

Austria: (0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20
Germany: +49 2304 932-5 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478
Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00
Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303
Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500